

Antibiotics used for selection pressure in culture media

Antibiotic	Stock	Final concentration
apramycin	100 mg/mL (in dH ₂ O)	50 µg /mL
hygromycin	50 mg/mL (in dH ₂ O)	100 µg/mL
kanamycin	50 mg/mL (in dH ₂ O)	50 µg /mL
chloramphenicol	25 mg/mL (in 100% EtOH)	25 µg /mL
carbenicillin	100 mg/mL (in dH ₂ O)	100 µg /mL
thiostrepton	50 mg/mL (in 100% DMSO)	50 µg /mL
nalidixic acid	25 mg/mL (in 0.1 M NaOH)	25 µg /mL